

A troubleshooting guide:
Real-time PCR experts
share their advice for
handling your top
experimental challenges

Real-Time PCR Tech Guide volume one

&A:

MX3005P™ System

Most Flexible

MX3000P® System

Most Affordable



Performance runs in the family.
Choose the personal QPCR system that's right for you.

Stratagene now offers two affordable, fully-featured quantitative PCR (QPCR) systems. The new five-color Mx3005P™ QPCR System includes expanded features to support a wider range of real-time QPCR applications, such as simultaneous five-target detection and alternative QPCR probe chemistries. The Mx3000P® QPCR System is still the most affordably priced four-color 96-well system available.

- A four- or five-color instrument, with user-selected filters
- Advanced optical system design for true multiplexing capability, and wider application support
- QPCR Software with enhanced data analysis and export functionality

Need More Information? Give Us A Call:

Stratagene USA and Canada

Order: (800) 424-5444 x3

Technical Services: (800) 894-1304 x2

Stratagene Europe

Order: 00800-7000-7000

Technical Services: 00800-7400-7400

Stratagene Japan K.K.

Order: 03-5159-2060

Technical Services: 03-5159-2070

www.stratagene.com

Mx3000P® is a registered trademark of Stratagene in the United States.
Mx3005P™ is a trademark of Stratagene in the United States.

Purchase of this product is accompanied by a license under the foreign counterparts of U.S. Patent Nos. 4,683,195, 4,683,202 and 4,965,188 covering the Polymerase Chain Reaction ("PCR") process, where such process is covered by patents. This instrument is an Authorized ThermalCycler for use with applications licenses available from Applied Biosystems. Its use with Authorized Reagents also provides a limited PCR license in accordance with the label rights accompanying such reagents.






Table of contents

Letter from the Editor 5

Index of Experts 5

Q1: How do you select and validate **normalization factors**? 7

Q2: What are your tips for improving **assay design**? 11

Q3: How do you choose or develop **standard protocols**? 13

Q4: What do you do to achieve better **primer specificity**? 15

Q5: How do you **enable amplification** from damaged or impure samples? 16

List of Resources 18

One enhancement

makes all the difference

ABgene's new **Absolute™ MAX RT kits** feature a unique enhancer, which eliminates double stranded DNA contamination and alleviates the need for DNase treatment. Consistently detecting less than 10 copies of RNA template, kits are optimized for compatibility with most thermocyclers, and are the largest off the shelf range on the market. For more information, visit Abgene.com, or call the sales team on 1 800 445 2812, quoting ref: GTRT for introductory discount. Alternatively, email sales@us-abgene.com or sales@abgene.com (rest of world).

Cat.No.	Product Description	Quantity	Cat.No.	Product Description	Quantity
AB-1369/A	Reverse-it™ MAX 1st Strand Synthesis Kit	40 rxn	AB-1332/A	Absolute™ MAX QRT-PCR Mix	100x50µl rxn

Real-time, Q-PCR Multiplexing

Using CAL Fluor[®], Quasar[®] and Black Hole Quencher[®] Dyes

Compatible with: ABI Prism[®] ■ Bio-Rad iCycler[®]
Cepheid SmartCycler[®] ■ Corbett Rotor-Gene[™]
Roche LightCycler[®] ■ Stratagene MX[®]

Multiplexing real-time, Quantitative PCR (rt, Q-PCR) assays using unoptimized reporter dye sets places unusual performance and data analysis demands on today's real-time, Q-PCR thermocyclers.

Biosearch has developed two NEW series of performance-optimized dyes for multiplex real-time, Q-PCR – the CAL Fluor and Quasar Dyes. Having emission maxima ranging from 522 nm to 670 nm, these new dyes are high performing, low cost alternative dyes for seven fluorophores commonly used on real-time thermocyclers:

TET[™], JOE[™], VIC[™], HEX[™],
Texas Red[®], Cy3[™] and Cy5[™].



Coupled with our BHQ[®] dark quenchers you now can realize the full potential of your Q-PCR thermocycler!

For a FREE Guide to PCR Multiplexing on your thermocycler contact us or visit
www.multiplexqpcr.com



81 Digital Drive, Novato, CA 94949-5750 USA
PHONE: +1.415.883.8400 FAX: +1.415.883.8488
www.multiplexqpcr.com
info@biosearchtech.com www.biosearchtech.com

A MATCH MADE IN HEAVEN



©2005 Biosearch Technologies, Inc. Black Hole Quencher, BHQ, CAL Fluor, and Quasar are registered trademarks of Biosearch Technologies, Inc., Novato, CA. TET, JOE, VIC and HEX are trademarks of Applied Biosystems, Texas Red is a registered trademark of Molecular Probes, "Cy" and "Cy Dye" are trademarks of Amersham Biosciences. "LightCycler" is a registered trademark of Roche Molecular Systems, Inc., Alameda, CA. "Prism" is a registered trademark of Applied Biosystems, "Cycler" is a registered trademark of Bio-Rad Laboratories, "SmartCycler" is a registered trademark of Cepheid, "Rotor-Gene" is a trademark of Corbett Research, "MX 3000P" and "MX4000" are trademarks of Stratagene. PCR is a proprietary technology covered by several US patents including US Patent Nos. 4,683,195, 4,683,202 and

4,965,188, and by issued and pending counterparts outside the U.S. These patents are owned by Roche Molecular Systems Inc. and have been sub-licensed by PE Corporation in certain fields. Depending on your specific application you may need a license from Roche or PE to practice PCR. Additional information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Roche Molecular Systems, 1145 Atlantic Avenue, Alameda, CA 94501 or Applied Biosystems business group of the Applied Corporation, 850 Lincoln Centre Drive, Foster City, CA 94404. In addition, the 5' nuclease assay and other homogeneous amplification methods used in connection with the PCR process may be covered by U.S. Patents 5,210,015 and 5,487,972, owned by Roche Molecular Systems, Inc. and by U.S. Patent 5,838,848, owned by The Perkin-Elmer Corporation. The purchase of Biosearch Technologies products does not, either expressly or by implication, provide a license to use this or other patented technology. Licensing information can be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, CA 94404 or the Licensing Department at Roche Molecular Systems Inc., 1145 Atlantic Avenue, Alameda, CA 94501.

Letter from the editor



Welcome to the inaugural issue of *Genome Technology's* reference guide series. This slim volume contains our first installment of answers to quantitative PCR problems.

Despite its youth, quantitative PCR has become ubiquitous in any lab dealing in nucleic acid quantification. The technology is conceptually simple, easy to use, and widely applicable. But these virtues belie a host of challenges, from improving assay design to working with degraded samples. These topics and more are addressed by our international panel of experts. Their field-tested tips are

gathered here to help you troubleshoot your way through any experimental scenario involving real-time PCR. The guide itself is user-friendly, sleek, and complete with a resource list based on expert advice highlighted in the text. In short, it's made to keep on the bench, preferably within easy reach of your lab's qPCR work station.

This guide owes its existence to the generosity of our contributors, whose time, expertise, and counsel have transformed our initial questions into an accessible technical manual. With special thanks to them, we hope you enjoy the following pages.

— *Jennifer Crebs*

Index of experts



Pamela "Scottie" Adams
Director, Molecular Biology Core Facility
Trudeau Institute



Mikael Kubista
Professor, Head of R&D
TATAA Biocenter



Claus Lindbjerg Andersen
Associate Professor
Molecular Diagnostic Laboratory
Department of Clinical Biochemistry
Aarhus University Hospital



Shawn Levy
Assistant Professor
Vanderbilt University



Vladimir Benes
Head, Genomics Core Facility, Scientific
Core Facilities, Services & Technologies Unit
European Molecular Biology Laboratory



Jo Vandesompele
Center for Medical Genetics
Ghent University Hospital



Stephen Bustin
Professor of Molecular Science
Queen Mary University of London




Michael Zianni
Manager/Research Associate
Plant-Microbe Genomics Facility
Ohio State University

*Tired of doing
gene expression analysis
one gene
at a time?*

GenomeLab™ GeXP will not only boost your performance, but your spirits, too. With high-throughput automation, flexible genetic analysis and eXpress Profiling™, you can achieve quantitative, cost-effective gene expression analysis with up to 30 genes in a single reaction. See why things are looking up at www.beckmancoulter.com/GeXP.





How do you select and validate normalization factors?

We test a panel of normalization genes — GAPDH, β -actin, β 2-microglobulin, GUS, and HPRT — to convince ourselves that expression is stable in the tissues and conditions we are using. GAPDH appears to be very stable in the tissues that we examine during various infections in the mouse. GAPDH is also useful because it can be used as an indication of genomic DNA contamination when used as a no transcription (-RT) control.

— *Pamela Adams*

Selection of candidates: I choose a set of candidate genes, typically five to 10, which I expect not to show expression variation under the conditions in the given experiment. Optimally the genes are selected based on some form of evidence (e.g. microarray data), but this is not a prerequisite (it will most likely increase the chance of identifying a good stable expressed gene among the candidates).

Validation of expression stability: I then purchase or design real-time RT-PCR assays to these genes and test these on samples identical to, or at least very similar to, the samples in the given experiment. For clinical samples I have a minimum of eight samples per sample group. I recommend that the generated C_t values are transferred to linear values using standard curves. If standard curves are not applied, caution should be taken to verify equal amplification efficiencies of all assays; alternatively, measures should be taken to correct for differences in amplification efficiencies.

I then estimate the intra- and inter-group expression variation of the individual assays using

the free-ware **NormFinder**. Finally, I select the candidate with the least expression variation as my normalization gene.

Alternatively, one can use a combination of two to three of the candidates to form a normalization factor — make sure that the candidates included have opposite directed inter-group variation. Otherwise the normalization factor is no better than using the candidates individually!

— *Claus Andersen*

Microarray data represent a treasure trove of reference genes for organisms where such data are available. Papers out there talk about genes expressed in normal tissues, and for some people that it is a very useful clue, but there are people who are working with organisms that are not so well characterized by microarrays. For these people, there is only testing, going through the literature, and more testing. Testing requires first selecting panels of reference genes: run them in a few situations, see if they are stable. This can be cumbersome, and there is the trouble that people are prone to compromise C_t difference, although they really shouldn't. This is a real issue that will need debate.

Some people share [Jo Vandesompele's] view regarding normalization factors, but it must be remembered that there is a difference between the ideal world and the world of daily research. There will be no such thing as a general set of reference genes. People need to be open minded and look to other biochemical pathways that are unrelated to the process they are *(continued on page 9)*


Real Time PCR --That Works!



SYBR® Premix *Ex Taq*™ (Perfect Real Time)

SYBR® Premix *Ex Taq*™ (Perfect Real Time) is a real time PCR premix from Takara which is unsurpassed for high-speed, high-specificity application results. This 2X concentration premix combines the high sensitivity and performance offered by Takara's *Ex Taq*™ Hot Start DNA Polymerase, SYBR® Green I, plus a newly formulated R-PCR buffer which provides superior specificity and increased amplification efficiency. Your result is exceptional real time PCR results quickly and easily.

- **Versatility:** use on many real-time PCR instruments such as: Smart Cycler®, LightCycler®, Applied Biosystems qPCR systems, MJ Opticon®, iCycler® and RotorGene™
- **High Sensitivity:** as few as 100 copies using SYBR® Green I detection.
- **Wide Dynamic Range:** this product possesses a dynamic range of 7-8 orders (λ DNA template).
- **High Speed Capability:** works with high speed qPCR instruments.
- **Convenient:** supplied in an easy-to-use format with the SYBR® Green I dye.



Performance of SYBR® Premix *Ex Taq*™ (Perfect Real Time) (left panel) vs. Roche's Fast Start DNA Master SYBR® Green I (right panel) using a Roche LightCycler®. Takara's SYBR® Premix *Ex Taq*™ demonstrates an excellent melting curve, particularly when low copy number templates are amplified. For more data demonstrating SYBR® Premix *Ex Taq*™ performance on other instruments, please see our website at www.takaramirusbio.com.



Takara
Takara Mirus Bio

www.takaramirusbio.com
888-251-6618

Order Your FREE Sample Today!

Ex Taq™ is a trademark of Takara Bio Inc. SYBR® is a registered trademark of Molecular Probes, Inc. Lightcycler® is a registered trademark of a member of the Roche group. RotorGene™ is a trademark of Corbett Research. Smart Cycler® is a registered trademark of Cepheid. iCycler and MJ Opticon are trademarks of BioRad laboratories Inc. Takara PCR Related Products are sold under a licensing arrangement with Roche Molecular Systems and F. Hoffman La Roche Ltd. and Applied Biosystems. Takara Bio's Hot-Start PCR-Related products are licensed under U.S. Patent 5,338,671 and 5,567,287 and corresponding patents in other countries.

How do you select and validate normalization factors?

(continued from p.7)

studying. Unfortunately, we are experiencing sort of a legacy issue with GAPD, actin, 18S genes, etc., because these genes were used at the beginning and now everyone is set to use them even when there is a lot of evidence out there showing that those genes can also be heavily regulated. Microarrays provide much richer data in this regard.

— *Vladimir Benes*

We have validated the use of total RNA in our particular experimental system (the colon) and express copy numbers per microgram total RNA. We have also developed a genomic DNA-based normalization method for RNA extracted from tiny samples, such as laser capture microdissected tissue. This uses Larry Wangh's PurAmp to lyse the tissue and perform the RT-PCR assay without having to purify the RNA any further (Hartshorn *et al.*, 2003).

— *Stephen Bustin*

We always test a large number of candidate reference genes, for example, using our endogenous control gene panel on representative samples, and then we select those with the least internal variation as references for the particular system. The latter selection can be done using, for example, **geNorm** or **GenEx**.

— *Mikael Kubista*

We either use several factors or genes for normalization to maximize the likelihood that one or more of the factors are appropriate, or experiments are run to empirically determine which available factors will be useful for the experiments at hand. In situations where RNA quality may be a challenge, we use the latter approach. This allows us to test several samples and the normalization factors at the same time.

— *Shawn Levy*

“People are prone to compromise C_t difference, although they really shouldn't. This is a real issue that will need debate.” — *Vladimir Benes*


According to our own established **geNorm** concept, we initially evaluate the expression stability of 10 candidate reference genes belonging to different functional and abundance classes. Using the **geNorm** software, we select the

optimal number (at least three) of most stably expressed reference genes and use the geometric mean of these reference genes as an accurate normalization factor.

— *Jo Vandesompele*

The customers select the standardization genes and then do a thorough optimization protocol and test to determine whether the C_t values are changing for a given gene across the samples. I also encourage them to use PCR products or cloned genes to check efficiency.

— *Michael Zianni*



Assays That Play Well Together

MGB
eclipse
PROBE SYSTEMS
For Real Time PCR

Anyone can design a Real Time PCR assay that 'plays alone.' The real challenge is developing multiplex assays where all components play well together. Bi-plexing a reference gene as an internal control with a target assay allows correction for sample to sample variation. Bi-plexing cannot be an afterthought however - it must be designed into the individual assays. Nanogen's MGB Eclipse™ By Design Service offers:

- Custom development of bi-plex assays to your specific target sequences
- Assays for targets not possible with other chemistries
- Elimination of lengthy and costly optimization steps
- Universal reaction conditions over a broad range of expression levels

And as any child knows, playing well together also means sharing.

Nanogen provides the sequence of your primers and probes to you. After all, they are your sequences.

NANAGEN

www.nanogen.com — 1-877-nanogen

Clinical Assays


Gene Expression

SNP Genotyping

Detection Assays

MGB Eclipse

This product is covered by one or more of the following US patents: 5,801,155; 6,492,346; 6,727,356; 6,030,787; 5,723,591; 5,876,930 and pending applications. For Research use only. The Polymerase Chain Reaction (PCR) is covered by US patents which are owned by F. Hoffmann-La Roche Ltd.



What are your tips for improving assay design?

I like to design three sets of primers for each assay, test them using a SYBR green assay, and then use the best set with a probe in a Taqman assay. If you design primers to span an intron to eliminate genomic DNA amplification, be sure to check the size of the intron. Some introns are small and if the amplicon amplified from the genomic DNA is only 100 bp bigger than an amplicon from mRNA, it could have an effect on the accuracy of your assay. Search the literature carefully to be aware of alternate splicing scenarios and what effect they might have on the assay.

— *Pamela Adams*

If possible, change the exon-intron boundary.

Redesign primers — be careful to avoid polymorphic and unspecific primer targets.

Instead of pairing primers based on melting temperature, primers should be paired according to priming efficiency (can be estimated using **OLIGO Primer Analysis Software**).

— *Claus Andersen*

I would say not to be greedy. These days, the price of oligos has dropped to point where it is not an issue to have two to three primer pairs in one or even different regions of transcribed sequence, which can be tested. One doesn't have to run a whole set of 96 reactions, but maybe a couple of reactions using each primer pair. In this way, you can proceed rather rapidly. The advantage of this approach is that if bad goes to worst, primers can be combined within primer pairs. Basically, we try to get nested primers if possible, so if one pair is not working, we can take the other primer

reversed from the amplicon to see if there has been an improvement.

There is this one thing which I have used recently by Premier Biosoft: they have software packages for microarray data, and **Beacon Designer** is the program for looking at secondary structures in the template when designing primers and probes. This software is a very nice tool, and there are also some Web pages where one can do similar things with the way sequences can fold.

— *Vladimir Benes*

Strive for maximum amplification efficiency; keep amplicons as short as possible; use specific primers where possible, otherwise oligo-dT; run standard curves with each assay; always run duplicate, no template controls, one set with tubes sealed prior to dispensing template, the other at the end of the dispensing. Always study amplification plots.

— *Stephen Bustin*

Use Evolution Operation (EVOP) optimization. If the assay is still not good enough, then change primers. For important assays we always test several primer combinations.

— *Mikael Kubista*

Have a good understanding of what you are designing and what the important considerations are. Like most things, the more you understand the process and what you are designing, the better the results will be. Also, practice makes perfect, so the more design and *(continued on page 17)*




- Easily programmable temperature gradient



- Optional motorized heated lid



- GLP-compliant environment



- Three different block formats



Licensed
and authorized
for PCR

Get your lab up to speed!

A new tempo that defies speed limits—rapidly evolving PCR technology meets its match with the Mastercycler® ep line of high-speed thermal cyclers.

Mastercycler ep is one of the fastest instruments in its class. A combination of high ramping speeds (up to 6 °C/s) and outstanding temperature regulation accuracy ensure extremely precise and fast PCR performance.

The new Mastercycler ep system:


- Extremely fast heating and cooling
- Standard PCR run times of ~30 min
- Intuitive graphic programming
- Superior temperature control technology
- No restrictions—use kits, reagents and consumables of your choice
- Remarkably quiet operation

For more information visit www.mastercycler.com

eppendorf
In touch with life

www.eppendorf.com • Email: info@eppendorf.com • Application hotline: 516-515-2258

In the U.S.: Eppendorf North America, Inc. 800-645-3050 • In Canada: Brinkmann Instruments (Canada) Ltd. 800-263-8715



How do you choose or develop standard protocols?

The protocol is dictated by the science. What question do you want to answer? I would search the literature for situations in which others have asked similar questions and see what protocols have been successful. Stephen Bustin's book, **A-Z of Quantitative PCR**, is a wonderful resource and **Michael Pffaf'l's webpage** contains a wealth of information.

— *Pamela Adams*

Roughly, the qPCR protocol can be divided into three parts: isolation of specimen and total RNA from it, a reverse transcription step, and a quantitative PCR step. Nowadays, there is no doubt that the last step is perhaps the most straightforward from the three thanks to numerous high-quality kits available on the market. It is not to say that one shouldn't work carefully, but not so many things can go wrong here, provided that one is approaching it with the right attitude. It is the first two steps which can considerably influence the outcome of the experiment.

Basically, the protocol should be very robust in terms of performance. Although ideally one would like to have one protocol that fits all, there are groups of protocols that may be more applicable. That is, if you know that some amplicons of some RNAs are folding close to 3', in which case you cannot have amplicon designed there, then it is necessary to move to 5' end and adjust your strategy of priming the RT reaction. Ultimately, regardless whatever protocol has been chosen, the most important factor is consistency within the cohort of

analyzed samples.

Preferentially, for any type of priming, the RT polymerase should be as thermostable and processive as possible. Mostly we are using SuperScript III and testing transcripts by Roche to see reproducibility over many runs. It is important to try to work with material that is as close to material one will get in a real experiment as possible. It is also better to use cDNA, not the plasmid, as the latter is too optimal. It is important to find the right measure, for too much of material can be inhibiting, for example.

We use a two-step method because the RT product can be archived; Stephen Bustin uses transcript-specific primers, but the product of this reaction cannot be used for assaying any other transcript, unfortunately.

Basically, people tend to think that qPCR is [an] easy method — it looks so deceptively simple — but it really needs to be recognized that there are risks and caveats involved.

— *Vladimir Benes*

We have developed our own protocols:

(a) In our hands, the most reliable quantification of RNA is achieved using target-specific primers for the RT step. We empirically choose an optimal primer pair from four potential candidates for each target gene.

(b) We use Mfold to determine the best location for reverse (RT) primers and usually carry out the RT step at 60°C for three minutes using *Tth* polymerase. We also use Mfold to determine the optimal position for the forward primer. *(continued on page 17)*



Roche Applied Science
LightCycler® 480 Real-Time PCR System

Rapid by nature,
accurate by design



▲ LightCycler® 480 Thermoblock for 96 or 384 wells, easily exchanged by users within minutes.

For years, Roche has provided real-time automated PCR solutions you can count on. Now, you can obtain the proven performance and benefits of the original LightCycler® System in a 96- or 384-well instrument platform for high-throughput applications – the new **LightCycler® 480 Real-Time PCR System**.

Speed – Save time without sacrificing the quality of your results – precise, high-speed temperature changes maximize specificity and yield.

Accuracy – Benefit from our novel thermal block and data-capture technologies to minimize edge-effects for outstanding accuracy and precision.

Versatility – Combine 5 excitation and 6 detection channels, multiple probe formats, proven analysis software, and true master mix reagents to meet your specific application needs.

Compatibility – Take advantage of the instrument's automation and LIMS capabilities to interface with your current systems and future workflows.

Visit www.roche-applied-science.com/lightcycler480 for more information.

In Development. Planned introduction: September, 2005.
For general laboratory use.
Not for use in diagnostic procedures.

This LightCycler® 480 Real-Time PCR System is licensed under U.S. Patent No. 6,814,934 and corresponding claims in its non-U.S. counterparts and under one or more of U.S. Patents Nos. 5,038,852, 5,656,493, 5,333,675, or corresponding claims in their non-U.S. counterparts, for use in life science research. It is also an Authorized Thermal Cycler. Purchase and use of this LightCycler® Instrument, in conjunction with Authorized Reagents, provides a limited license for use of the PCR process in life science research. No rights are conveyed expressly, by implication or by estoppel under any other patent claims or for any other application.

LIGHTCYCLER is a trademark of Roche.


The technology used for the LightCycler® System is licensed from Idaho Technology, Inc., Salt Lake City, UT, USA.

© 2005 Roche Diagnostics GmbH. All rights reserved.



Diagnostics

Roche Diagnostics GmbH
Roche Applied Science
68298 Mannheim
Germany



What do you do to achieve better primer specificity?

The best way to achieve primer specificity is by careful assay design. I prefer the extra assurance of specificity provided by the Taqman assay.

— *Pamela Adams*

Although it helps to Blast the primers to see if there is any secondary binding within transcriptome tested, people are not always working with organisms where sequence information is available. I think the ultimate test of your primers is a standard RT-PCR and checking its outcome on the high-percentage agarose gel. If the result is satisfying, we're running a small sample set on [the] qPCR instrument with SYBR Green and look at the dissociation curve. If the primer pair passes this test, then it's suitable for work.

— *Vladimir Benes*

Take nucleic acid secondary structure into account, perform the RT step at 60°C, reduce primer concentration as much as possible, and if all else fails, reverse transcribe from mRNA.

— *Stephen Bustin*

We don't experience primer specificity to be a problem, but poorly designed primers may give rise to primer-dimer products. This is minimized by proper primer design including avoiding any 3' complementarity, using low primer concentrations, hot-start, and short annealing times.

— *Mikael Kubista*

Primer specificity is best achieved through proper design of the assay. Trying to optimize conditions to get a sub-optimal primer set to perform at high efficiency is typically more problematic than doing a good design or using a validated design to start with.

— *Shawn Levy*


Limit GC content of the 3' end of the primer (maximum 40 percent GC content for last five bases). After primer design, we routinely do a Blast specificity search using the 'search for short, nearly exact matches' option at **NCBI's Blast portal** and using both forward and reverse primer sequence separated by three N's as one input sequence.

For more difficult designs where there is less freedom to choose the sequence to amplify (e.g. SNP genotyping or splice variant quantification), we sometimes add agents to enhance specificity, such as DMSO, or increase annealing temperature. We also keep the MgCl₂ and primer concentrations limited (max 3.5 mM and 250 nM for SYBR Green I assays, respectively), as higher concentrations allow more mispriming.

— *Jo Vandesompele*

Most customers use SYBR Green and have limited budgets, so they use thermal cycle changes instead of ordering different primers. Most customers (1) raise the annealing temperature, (2) collect data at a temperature intermediate to extension and melting, and/or (3) use five to 10 seconds for annealing.

— *Michael Zianni*



How do you enable amplification from damaged or impure samples?

If your amplicon is really short, then you still can use the sample and get decent curves, even if the mRNA has been heavily degraded. However, cDNA can be obtained only by priming the RT reaction with random primers. It's always good to check using the **Agilent Bioanalyzer** to see the profile of the sample's total RNA.

One thing that has been introduced in the lab is an assay assessing the scale of degradation. In regard to this, I recommend the work of Tania Nolan and Stephen Bustin. They have introduced two assays that are useful. One is the 3':5' strategy, that is checking quality of the 3' and 5' ends of selected transcript by qPCR and look to see if there is any C_t curve difference. If the mRNA is intact, then you should get the same C_t values for both ends. Nolan and Bustin also invented a method to assess presence of inhibitors; nevertheless, people try to push their luck here.

So, if a sample is very precious, with no chance of recovery and such that you cannot apply any further purification steps, then the length of the amplicon is very important. If the issue is degradation, you cannot do oligo-dT priming, but are priming with random primers, it is important that it is consistent with all samples in a panel. In any case, there is always hope in these situations (of course, if it cannot be done, then sorry).

— *Vladimir Benes*

It is important to realize that results from impure/damaged RNA must be treated with caution. We screen our RNA using the **Agilent Bioanalyzer**,

and assess the 3':5' ratios of reference genes, and screen for inhibitors present in the sample.

— *Stephen Bustin*

When amplifying damaged samples we believe one should use specific RT primers, which are the same as one set of the PCR primers, to minimize variation in the RT step due to RNA degradation.

For complex samples it is important to account for matrix effects, which can be done either by *in situ* calibration (standard dilution) or standard addition (Stahlberg *et al.*, 2003).

— *Mikael Kubista*

RNA isolation procedures that are optimized for the conditions of the sample (frozen or FFPE, tissue or cell lines, etc.) make a tremendous difference in the quality of RNA. When the quality of the RNA cannot be made better, ensuring that all of the RNA that will be used in differential comparisons has the same overall quality is very important to ensure consistent and comparable results. If impurities are present, cleaning up the RNA samples prior to cDNA synthesis or cleaning up the cDNA material can be very helpful.

— *Shawn Levy*

We have experienced that it is fairly easy to obtain beautiful amplification plots for degraded samples. The important question, however, is how reliable data are coming from such samples. In a recent study, we demonstrated that expression patterns might vary in relation to the degradation status of

the samples (Perez-Novo *et al.*, 2005). In the same paper, we recommended not to analyze intact and degraded samples together.

With respect to working with impure samples, we try to address this issue by performing proper nucleic acid extraction, followed by extensive quality

control, such as inspecting the UV absorption spectrum of the nucleic acids and performing capillary gel electrophoresis to assess the length of the ribosomal RNA molecules.

— *Jo Vandesompele*

Q2: What are your tips for improving assay design?

(continued from p.11)

testing, the better you get. Of course, with the number of vendors selling pre-validated assays, the easiest way to improve your assay design is to buy one rather than build one.

— *Shawn Levy*

Perform extensive *in silico* assay evaluation, using Blast specificity analysis of the primer pair, Mfold secondary structure prediction of the amplicon sequence (significant structures reduce amplification efficiency), search for SNPs present in the primer annealing region, and transcript variant analysis (alternative exons, etc.).

We are currently implementing a Web-based and freely available *in silico* assay evaluation analysis pipeline for user-designed assays. This module will be part of the next update of our real-time PCR primer and probe database (**RTPrimerDB**) and aims at streamlining assay quality prediction prior to experimental evaluation, based on above cited established *in silico* tools, such as Blast specificity search, Mfold secondary structure analysis, presence of SNPs or plain sequence errors, and graphical visualization of the aligned primer sequences on the target gene.

— *Jo Vandesompele*

Do a thorough optimization at the start focusing on annealing temperature, template concentration, primer concentration, and utilizing the best bench technique as possible while setting up the reactions. I make sure customers realize there are no shortcuts or magic bullets, and thorough optimization is critical.

— *Michael Zianni*

Q3: How do you choose or develop standard protocols?

(continued from p.13)

(c) A standard curve is developed for each target. We aim to achieve an amplification efficiency of more than 95 percent. If that is not possible, we aim to achieve consistent amplification efficiencies between runs, monitored using a standard curve as well as a positive control sample.

— *Stephen Bustin*

When using standard protocols, it is important that all assays are optimized under standard conditions. If this is not the case, matrix effects from biological samples will complicate interpretation seriously.

— *Mikael Kubista*

We start with recommended or established protocols and then make modifications as dictated by experimental needs or try optimizations to improve performance or efficiency. For the most part, we use established and well used protocols — usually from the manufacturer. If it isn't broken, don't fix it.

— *Shawn Levy*

Standard protocols or operating procedures are very important in order to achieve consistent and accurate results. As our lab was one of the early adopters of qPCR technology, we developed our own experimental protocols, from nucleic acid quality control, over DNase treatment, cDNA synthesis, PCR reaction setup, and data analysis. It is important to rigorously control the quality at each step and use both positive and negative controls to evaluate the reliability of the results.

— *Jo Vandesompele*

List of resources

In the preceding pages, our experts referred to a number of resources. We've compiled their recommendations in the following list. Whether you're a qPCR ace or apprentice, these Web tools and publications are sure to be helpful.

Publications

Andersen CL, Ledet-Jensen J, Ørntoft T (2004). Normalization of real-time quantitative RT-PCR data: a model based variance estimation approach to identify genes suited for normalization, applied to bladder and colon cancer data sets. *Cancer Research*, 64: 5245-5250.

Bustin SA (2004). *A-Z of Quantitative PCR*. International University Line: La Jolla.

Bustin SA, Benes V, Nolan T, Pfaffl MW (2005). Quantitative real-time RT-PCR - a perspective. *J Mol Endocrinol*, 34(3): 597-601.

Bustin SA, Nolan T (2004). Pitfalls of quantitative reverse transcription polymerase chain reaction. *J Biomol Tech*, 15: 155-166.

Hartshorn C, Anshelevich A, Wangh LJ (2005). Rapid, single-tube method for quantitative preparation and analysis of RNA and DNA in samples as small as one cell. *BMC Biotechnol*, 5(1): 2.

Perez-Novó CA, Claeys C, Speleman F, Van Cauwenberge P, Bachert C, Vandesompele J (2005). Impact of RNA quality on reference gene expression stability. *Biotechniques*, 39(1): 52, 54, 56.

Stahlberg A, Aman P, Ridell B, Mostad P, Kubista M (2003). Quantitative real-time PCR method for detection of B-lymphocyte monoclonality by comparison of kappa and lambda immunoglobulin light chain expression. *Clin Chem*, 49(1): 51-59.

Vandesompele J, De Preter K, Pattyn F, Poppe B, Van Roy N, De Paepe A, Speleman F (2002). Accurate normalization of real-time quantitative RT-PCR data by geometric averaging of multiple internal control genes. *Genome Biol*, 3(7): RESEARCH0034.

Websites

Agilent Bioanalyzer: <http://www.chem.agilent.com/>

GenEx: <http://www.multid.se/GenEx/genex.htm/>

Gene Quantification web page (edited by Michael W. Pfaffl): <http://www.gene-quantification.info/>

GenNorm: <http://medgen.ugent.be/~jvdesomp/genorm/>

NCBI's Basic Local Alignment Search Tool: <http://www.ncbi.nlm.nih.gov/blast/>

NormFinder: <http://www.mdl.dk/publicationsnormfinder.htm/>

OLIGO Primer Analysis Software: <http://www.olygo.net/>


Premier Biosoft Beacon Designer: http://www.premierbiosoft.com/molecular_beacons/

Real Time PCR Primer and Probe Database (RTPrimerDB): <http://medgen.ugent.be/rtprimerdb/>

TATAA Biocenter: <http://www.tataa.com>

Optimized, Flexible, Sensitive...

We Have the Right QPCR Reagent for You.



Sigma's broad range of QPCR reagents provide the flexibility needed to meet your QPCR needs. From probe based to SYBR® Green dye based QPCR, our products are compatible with all instruments and optimized to help you achieve successful results.

- Optimized ReadyMixes™ let you choose the detection method.
- Compatible with tube or plate based instruments.
- Greater sensitivity, efficiency, and linear response.

Probe based Quantitative PCR	
Product Number	Product Name
D6442	JumpStart™ Taq ReadyMix for High Throughput Quantitative PCR (blended with Reference dye)
D7440	JumpStart Taq ReadyMix for Quantitative PCR
D9191	JumpStart Taq ReadyMix with dUTP for Quantitative PCR
SYBR Green based Quantitative PCR	
Product Number	Product Name
S9194	SYBR Green JumpStart Taq ReadyMix for High Throughput Quantitative PCR (blended with Reference dye)
S4438	SYBR Green JumpStart Taq ReadyMix for Quantitative PCR
S5193	SYBR Green JumpStart Taq ReadyMix without Magnesium Chloride
Capillary based Quantitative PCR	
Product Number	Product Name
S1816	SYBR Green JumpStart Taq ReadyMix for Quantitative PCR, Capillary Formulation
D9191	JumpStart Taq ReadyMix with dUTP for Quantitative PCR
Quantitative RT-PCR	
Product Number	Product Name
QR0100	SYBR Green Quantitative RT-PCR Kit
QR0200	Quantitative RT-PCR ReadyMix for probe based applications

For more information and to evaluate Sigma's Quantitative PCR products, contact us at: sigma-aldrich.com/qpcr1 or call 1-800-325-3010

SYBR® Green nucleic acid stain is licensed from Molecular Probes, Inc. under U.S. Patent No. 5,436,134 for use in Real-Time PCR by end-users for research and development only. SYBR® is a registered trademark of Molecular Probes, Inc.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA






Pass the iQ Test

What do all these images have in common? The answer is 5. Introducing the iQ5 real-time PCR detection system and reagents, bringing true five-color detection capabilities to all of your real-time PCR applications.

Five reasons why the iQ5 is the ultimate system for value and flexibility:

- Powerful thermal gradient
- Optical detection system features user-changeable filter sets
- Gene expression analysis by relative quantity (ΔC_T) or normalized expression ($\Delta\Delta C_T$)
- Customizable software reports that offer comprehensive summaries of assay run conditions, data graphs and tables, and data analysis parameters
- Optimized multiplex reagents: increase throughput, control costs, and maximize your data output with the new iQ multiplex powermix



NEW

iQ5 real-time PCR detection system
and iQ™ multiplex powermix

For more information, visit us on the Web at
www.bio-rad.com/products/iQ5/ or
www.bio-rad.com/supermixes/