

Nikon Instruments Eclipse Ti-E Inverted Microscope coupled with Nipkow spinning disk confocal imager CARV II

Location

Room D / basement / 048 (l. 2472)

Institute of Physiology of the Czech Academy of Sciences (CAS), Areal Krč, Vídeňská 1083, Prague 4

Usage rules

See the document [General Usage Rules for Confocal Microscopes](#)

Specification of the confocal microscope Nikon Ti-E / Crest CARV II

1. **Microscope.** The inverted fluorescent microscope Nikon Ti-E with motorized stage Prior Proscan III enables one to conduct multiple ROI, Tile Scan and Mark&Find experiments. The microscope is coupled with Nipkow spinning disc imager CARV II that contains a unit for confocal and FRAP imaging. The whole system is placed on an active pneumatic anti-vibration table STable© Supertech.

2. Available objectives.

- Nikon CFI Plan Fluor 10X, DIC, 10x/0.3 NA, WD = 16 mm;
- Nikon CFI ADL 10X, Ph, 10x/0.25, WD = 6.2 mm;
- Nikon CFI Plan Fluor 20X, DIC, 20x/0.5 NA, WD = 2.1 mm;
- Nikon CFI S Fluor 40X Oil, DIC, 40x/1.3 NA, WD = 0.22 mm;
- Nikon CFI Plan Apo Lambda 60X Oil, DIC, 60x/1.4 NA, WD = 0.13 mm;
- Nikon CFI Plan Apo VC 100X Oil, DIC, 100x/1.4 NA, WD = 0.13 mm;

For more information, please visit [Nikon website](#)

3. **Illumination path.** The sample can be illuminated by a halogen lamp in transmission mode or by the mercury lamp X-Cite® 120PC Q in epi-fluorescent mode. The provided spectrum ranges from 300 nm to 700 nm and the illumination path is optimized for common dyes (DAPI, GFP, Cy3, Texas Red, Cy5).

4. **Detection path.** The emission path is equipped with band pass filters that are optimized for the above mentioned dyes. The image is acquired by a cooled

monochromatic 16bit CCD camera Hamamatsu Orca-R² with resolution 1344 (H) x 1024 (W), pixel size 6,45 µm x 6,45 µm, acquisition speed 16,2 fps and quantum efficiency over 70 %.

5. Box incubator. The system Okolab UNO-COMBINED-CONTROLLER controls the environmental parameters such as temperature (ambient temperature to 50° C), CO₂ (0 to 15 %) and relative humidity (up to 75 % for 37° C). The size of the box is 85,5 x 127,7 x 25,0 mm.

6. Software. The whole system is controlled through NIS-Elements AR that supports a tool for running complex tasks (module Jobs).

Useful websites

- [Reservation system](#) at the Institute of Physiology, Czech Academy of Sciences:
- [Specification of optical microscopes](#) at the Institute of Physiology, Czech Academy of Sciences:

- Detailed description of the microscope:
[Internal users](#)

[External users](#)

- Information on the equipment that has been included into the project [Czech-BioImaging](#)
- [Training protocol – internal users](#) (including the rules for entering the facilities with laser scanning confocal microscopes at the Institute of Physiology, Czech Academy of Sciences, [template PDF](#))
- [Training protocol – external users](#)
- Image processing in free software [NIS-Elements Viewer](#)
- Image processing and analysis with open source program Fiji (Fiji is Just ImageJ) that includes a number of useful plug-ins:
[Bioformats \(former LOCI Tools\)](#)
[SLIM Curve](#) or <https://slim-curve.github.io/>
[ImageJ world mailing list](#)
- Interactive dye spectra viewer:



Crest CARV II/Nikon Ti-E @ IPHYS CAS
Updated 6. 4. 2018

[ThermoFisher](#) (Life Technologies) interactive Spectra Viewer

[Leica FluoScout](#)

[BdBioSciences Spectrum Viewer](#)

[BioLegend SpectraAnalyzer](#)

[-Tables of fluorescent dyes spectra](#)