


Monitoring of the Energy Scale in the KATRIN Neutrino Experiment



Martin Slezák

Faculty of Mathematics and Physics, Charles University in Prague
Nuclear Physics Institute, Řež near Prague

Understanding Neutrinos
24th Indian-Summer School of Physics, 2012
Prague, Czech Republic



Outline

1. Introduction to topic


- neutrino mass
- measuring neutrino mass: KATRIN experiment
- high voltage in KATRIN - monitoring concept
- $^{83}\text{Rb}/^{83\text{m}}\text{Kr}$ source as nuclear & atomic standard

2. Measurement results

- investigation of new line-shape: Doniach – Sunjic function
- stability measurements

Neutrino mass


- non-zero neutrino mass **confirmed** by various **oscillation experiments** (Super-Kamiokande, SNO, Daya-Bay, ...) excluding other explanations
- neutrino mass hierarchy?
- **absolute mass scale?**
- motivation for absolute mass scale also from **cosmology** (hot dark matter?)



- **current accepted limit** on neutrino mass (from T_2 β -decay measurements in Mainz & Troitsk):

$$m(\bar{\nu}_e) < 2 \text{ eV}/c^2$$

- potential: improve sensitivity **down to 0.2 eV**




Measuring neutrino mass: β -decay

- kinematic determination of electron anti-neutrino mass from β -decay – model independent
- decay rate given by:


$$\frac{d\Gamma}{dE} = C p(E + m_e)(E_0 - E) \sqrt{(E_0 - E)^2 - m_{\nu_e}^2} F(Z + 1, E) \Theta(E_0 - E - m_{\nu_e})$$

$$m_{\nu_e}^2 = \sum_{i=1}^3 |U_{ei}|^2 m_i^2$$

- suitable β -emitter: **tritium**
 - $E_0 = 18.6$ keV, $T_{1/2} = 12.3$ y
 - super-allowed transition
 - simple final states
 - (**KATRIN** experiment)
- other option:
 - ^{187}Re , $E_0 = 2.47$ keV,
 $T_{1/2} = 43$ Gy
 (MARE experiment)




KARlsruhe TRitium Neutrino experiment




- β -electrons are transported along magnetic field lines from the T_2 source through the transport and pumping sections and spectrometers to the detector
- T_2 pumped out by TMP's and cryosorption pumps (flow reduction factor of 10^{14})
- most of electrons pre-filtered at pre-spectrometer
- low background at main spectrometer, extreme high vacuum (10^{-11} mbar)

Measuring neutrino mass: MAC-E filter for KATRIN



- principle: **adiabatic guidance of electrons** from T_2 source along magnetic field lines to a plane where their energy is analysed using an electrostatic potential = **MAC-E filter**
- adiabatic transport: $\mu = E_{\perp} / B = \text{const.}$
- B-field drops by 4 orders of magnitude from solenoid to the analysing plane: $E_{\perp} \rightarrow E_{\parallel}$
- only electrons with $E_{\parallel} > eU$ can pass the retardation potential \rightarrow **integral spectrum**
 - acquired by recording electron countrate for particular retarding voltages U in **discrete voltage steps**
- **superior energy resolution**: $\Delta E = E_{\perp, \text{max, start}} \cdot B_{\text{min}} / B_{\text{max}} = 0.93 \text{ eV}$

High voltage in KATRIN - monitoring concept




- retardation potential depends on **high voltage** applied to the spectrometer
- for precise analysis of electron energy by retardation potential the high voltage has to be known to **very high precision**
 - to achieve desired ν mass sensitivity HV has to be stable within ± 60 mV over 2-month KATRIN run \rightarrow for $U \sim 18.6$ kV: **± 3 ppm !**

Two methods of monitoring the HV:

- reduce HV by a **high-precision divider** and measure it directly with a **precise voltmeter**
- apply the same HV on **another MAC-E filter** (= **Monitor spectrometer**) and measure **time-stable conversion electron line** by varying potential on the source
 - any change in line position points to instability of the divider-voltmeter setup


KATRIN Monitor spectrometer



- former Mainz spectrometer, rebuilt in KIT Karlsruhe
- in practice: smaller version of the Main spectrometer
 - energy resolution 0.93 eV @ 18.6 keV
 - Earth & low magnetic field compensation system
 - ultra-high vacuum $\sim 10^{-9}$ mbar in the tank, $\sim 10^{-10}$ mbar in the source and detector chambers
 - 3D positioning of source and detector


Solid $^{83}\text{Rb}/^{83\text{m}}\text{Kr}$ conversion electron source

- best suited for monitoring of HV in KATRIN: **K-32 conversion electron line** from $^{83\text{m}}\text{Kr}$ – energy **800 eV** below tritium endpoint
- $T_{1/2}$ of $^{83\text{m}}\text{Kr}$ only 1.8 hours
→ **replenishment from ^{83}Rb** with $T_{1/2} = 86$ days
- method: to **capture $^{83\text{m}}\text{Kr}$ atoms** in Pt substrate before they decay
→ solid sources
- production: **implantation** of ^{83}Rb atoms into pure platinum substrate at ISOLDE mass separator (CERN)
- under investigation: implantator at HSKP Bonn




* D. Vénos et al., NIM A 560 (2006) 352

** M. Slezák et al., EPJ A 48 (2012) 12




Integral spectrum of the MAC-E filter: analysis



- the integral spectrum is analysed by means of **many-parameter fit** (search of chi-square minimum via Minuit routine under ROOT)
- **expected** conversion line shape: Lorentzian function
- **observed** conversion line shape: Lorentzian \otimes Gaussian & asymmetrical shape at lower energy side (Lor. \otimes Gaus. = Voigt shape)
- previously used description: **two Voigt lines** (motivation: ^{83}Rb atoms sitting in two different environments in the substrate)
- new approach: single **Doniach & Sunjic asymmetrical line shape** convoluted with Gaussian (motivation: interaction of conversion electrons with conduction electrons in the metal)
 - S. Doniach, M. Sunjic, J. Phys. C 3 (1970) 285


Line shapes – differential form examples



Voigt (symmetrical)

- no asymmetry was observed in Mainz for vacuum evaporated solid sources
→ has to be a solid-state effect
- one source was annealed up to 400 °C
→ no change in line shape


**Doniach-Sunjic (conv. with Gauss.)
(asymmetrical)**



asymmetrical part, for zero asymmetry
merges into Voigt


Line shapes – integral form examples (1)

One Voigt fit:




× does not work at all,
data are not described
well

Two Voigt fit:



One Doniach-Sunjic fit:




- works well, but only on **limited space**
- left peak tends to fluctuate
- **no. of free parameters: 6**

- ✓ works perfectly
- ← ✓ can describe more data
- ✓ **no. of free parameters: 5**

Line shapes – integral form examples (2)

- Doniach-Sunjic lineshape works with other conversion electron lines in the spectrum as well:


L_3 -32, source Pt10-1



sharp line (small lorentzian width)

broad line (large lorentzian width)


M_1 -9.4, source Pt10-1



- asymmetry comparable among different lines
 - was not the case with 2 Voigt fits


K-32 stability measurements (1)

- for the monitoring concept of KATRIN to work, the K-32 line energy has to be **stable in time** (max. drift ± 1.5 ppm per month)



K-32 stability measurements (2)

- possible explanation for large line drift: various experiments done at MoS introducing gases into spectrometer, discharges?



Background experiments at MoS

Venting & bake-out

Changes to HV system

Conclusions

- **absolute mass scale** of neutrinos unknown
→ measurement of the neutrino mass in the **KATRIN experiment** from T_2 β -decay
- **monitoring of high voltage** in KATRIN one of crucial points for reaching the designed sensitivity on the neutrino mass
- one method for monitoring: **another spectrometer** (Monitor spec.) connected to the same high voltage measuring **time-stable conversion electron line** from a solid $^{83}\text{Rb}/^{83\text{m}}\text{Kr}$ source
- Monitor spectrometer **successfully running** at KATRIN site at KIT
- K-32 energy stability **confirmed in two time intervals**, but between them the line position is **not as stable as required**
- possible reason: load to MoS due to **various experiments performed**
→ during KATRIN T_2 runs no such experiments will be done and conditions will be stable → **stable K-32 energy expected**
- **next major task**: demonstrate stability and monitoring capabilities in KATRIN mode (HV cable between spectrometers, varying potential on the source)

Thank you for attention!

