

Laboratory of Signal Transduction

Plasma membrane signalosomes, immunoreceptor signalling, mast cell activation, chemotaxis


Petr Dráber


petr.draber@img.cas.cz

www.img.cas.cz/research-groups/petr-draber

The main research interest of our laboratory has been traditionally focused on understanding molecular mechanisms governing signal transduction from the high affinity IgE receptor (FccRI) and other plasma membrane receptors to the cytoplasm. We continued in functional studies on the role of transmembrane adaptor proteins in regulation of such important antigen-mediated events as calcium response, degranulation and chemotaxis. Especially, we attempted to solve the problem we noticed in our previous studies that NTAL has a negative regulatory role on FccRI activation in mast cells isolated from NTAL knockout mice, but positive regulatory role in human mast cells with NTAL knockdown (through RNA interference). We found that the observed discrepancies did not reflect compensatory developmental alterations in mouse cells, as expected, but rather different roles of NTAL in FccRI-induced signalling pathways in human and mouse cells. To gain more insight into NTAL function, we also examined gene expression profiles in resting and antigen-activated mast cells with NTAL knockout and knockdown and corresponding controls and we identified several genes that were differentially expressed in NTAL deficient and wild-type cells. We also continued in studies on stromal interacting molecule 1 [STIM1] in mast cell signalling and found that STIM1 is required for formation of microtubule protrusions in antigen-activated cells and for chemotaxis towards antigen.

Several other molecules were analysed (Csk, PAG, STAT5, CD9, ORMDL3, FccRI) to get better insight into the signalling pathways involved in mast cell signalling. Significant effort was devoted to identifying new signalling molecules by high-throughput screening using RNA interference libraries.





Time (10 min)

Fig. 2. Dynamics of STIM1-Ora1 crosstalk. Interactions between STIM1 and Orai1 were determined by fluorescence resonance energy transfer [FRET]. FRET [in red] occurs between STIM1-EYFP and Orai1-ECFP in cells after activation with thapsigargin that releases calcium from endoplasmic reticulum. FRET is determined by flow cytometry.

- AS CR, M200520901 New components of immunoreceptor signaling in mast cells, 2009-2012, Pe. Dráber
- AV CR, KAN200520701 Nano-PCR ultrasensitive test for detection of specific proteins in body fluids, 2007–2011, Pe. Dráber
- GA CR, GA301/09/1826 Topography and function of Csk-binding proteins of the plasma membrane in mast cells, 2009-2013, Pe. Dráber
- GA CR, GAP302/10/1759 Function-structure relationships between transmembrane adaptor-based signalosomes in mast cells, 2010-2013, L. Dráberová
- Ministry of Industry and Trade of the Czech Republic, FR-TI3/067 Genetically modified polymerases and their use for amplification of unpurified DNA, 2011-2014, Pe. Dráber
- TA CR, TA01010436 New generations of DNA aptamers, 2011-2013, Pe. Dráber, L. Dráberová
- Ministry of Education, Youth and Sports of the Czech Republic, LD12073 LD-COST CZ Membrane signalosomes of mast cells and basophils targets for innovative therapies, 2012-2014, P. Dráber
- 1. Rashid A, Housden JEM, Helm BA, <u>Draber P</u>. Fc receptor- γ subunits with both polar or non-polar amino acids at position of T22 are capable of restoring surface expression of the high-affinity IgE receptor and degranulation in γ subunit-deficient rat basophilic leukemia cells. **Mol** Immunol 2012 53(3): 270-273.
- 2. Halova J, Draberova L, Draber P. Mast cell chemotaxis chemoattractants and signaling pathways. Front Immunol 2012 3: 119.
- 3. Hájková Z, <u>Bugajev V</u>, Dráberová E, Vinopal S, <u>Dráberová L</u>, Janáček J, <u>Dráber P</u>, Dráber P. STIM1-directed reorganization of microtubules in activated mast cells. **J Immunol 2011** 186(2): 913-923.
- 4. Horáková I., Polakovičová I., Shaik GM, Eitler J. Bugajev V. Dráberová L., Dráber P. 1,2-propanediol-trehalose míxture as a potent quantitative real-time PCR enhancer. BMC Biotechnol 2011 11: 41.
- 5. <u>Potůčková L, Franko F, Bambousková M, Dráber P</u>. Rapid and sensitive detection of cytokines using functionalized gold nanoparticle-based immuno-PCR, comparison with immuno-PCR and ELISA. **J Immunol Methods 2011** 371[1-2]: 38-47.



From the left:

Pavol Utekal, MSc / Research Assistant [since 2012] · Monika Bambousková, MSc / PhD Student · Tomáš Paulenda, MSc / PhD Student · Viktor Bugajev, MSc / PhD Student · Ivana Hálová, PhD / Research Fellow · Oleksii Redchenko, PhD / Research Fellow · Petr Dráber, DSc / Head of Laboratory · Lucie Potůčková, MSc / PhD Student · Lubica Dráberová, PhD / Research Fellow · Hana Mrázová / Technician, Secretary · Iva Polakovičová, MSc / PhD Student · Romana Budovičová, MD / Research Assistant

Not in the picture: Helena Dráberová, MSc / PhD Student · Lukáš Kocanda / Technician · Magda Tůmová, PhD / Research Fellow